Discoplan-TS
Instruction Manual

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
Struers A/S
Valhøjs Allé 176
DK-2610 Rødovre, Denmark
Telephone 36 70 35 00

Herewith declares that

Product Name:	Discoplan-TS
Type No:	372
Machine Type:	Grinding machine

is in conformity with the provisions of the following directives:


Safety of Machinery 89/392/EEC and 91/368/EEC and 93/44/EEC according to the following standard(s):
EN292-1:1991, EN292-2:1991, EN60204-1:1992.

EMC-Directive 89/336/EEC and 92/231/EEC according to the following standard(s):
EN50081-1:1991, EN50082-1:1992.

Low Voltage Directive 73/23/EEC and 93/68/EEC according to the following standard(s):
EN60204-1:1992.

Supplementary Information The equipment complies with the American standards:
UL508.

The above has been declared according to the global method, module A

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Discoplan-TS Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully instructed in the use of the apparatus according to the manual. The operator should be fully instructed in the use of cut-off wheels and grinding discs for Discoplan-TS.
2. Use only intact diamond cut-off wheels and grinding discs. If Struers cut-off wheels are not used, it should be noticed whether they contain materials which give off harmful vapors or dust requiring special ventilation. If Struers cut-off wheels are not used, the wheels must be approved for a spindle speed of minimum 1680 rpm.
3. Make sure to have suitable ventilation when cutting materials which give off harmful dust particles or vapors. See Safety Data Sheet for the materials in question.
4. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the cutting fluid in question.
5. Store the cutting fluid and use it as indicated in the instructions for the cutting fluid in question so that harmful bacterial growth and contact with the skin is avoided.
6. Always use safety goggles in connection with manual cutting. Use safety goggles when dressing the cut-off wheel or the grinding disc with a dresser stick. Use hearing protector when cutting.
7. Adjust the wheel guard according to the thickness of the specimen, so that it covers as much of the cut-off wheel as possible.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment. Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

Section 12: Other information

To be read carefully before use

1. The user must read the instructions carefully before using the product. The user must also read the instructions for the use of the product in the event of an emergency.

2. The user must read the instructions carefully before using the product. The user must also read the instructions for the use of the product in the event of an emergency.

3. The user must read the instructions carefully before using the product. The user must also read the instructions for the use of the product in the event of an emergency.

4. The user must read the instructions carefully before using the product. The user must also read the instructions for the use of the product in the event of an emergency.

5. The user must read the instructions carefully before using the product. The user must also read the instructions for the use of the product in the event of an emergency.

6. The user must read the instructions carefully before using the product. The user must also read the instructions for the use of the product in the event of an emergency.

7. The user must read the instructions carefully before using the product. The user must also read the instructions for the use of the product in the event of an emergency.

8. The user must read the instructions carefully before using the product. The user must also read the instructions for the use of the product in the event of an emergency.

Discoplan-TS Instruction Manual

IMPORTANT

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the front cover of this manual or on an attached label below. If in doubt consult the rating plate of the machine itself.

We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals

A Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals

A Struers Service Manual may only be used by a trained technician authorized by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 1996.

Struers A/S

Valhøjs Allé 176
DK-2610 Rødovre/Copenhagen
Denmark
☎ Telephone +45 36 70 35 00
☎ Fax +45 38 27 27 01

Table of Contents

Page

1.	Installation	5
2.	Operation	7
3.	Equipment, Consumables and Accessories	15
4.	Maintenance	19
5.	Technical Data	21

1. Installation

1.1 Packing List

- 1 Discoplan-TS bolted to the bottom of the transport case
- 1 Box spanner for bolts on the case, width 17mm
- 1 Key for mounting of diamond cut-off wheel, width 27mm
- 1 Hook-spanner for mounting of diamond cut-off wheel
- 1 Allen key width 2.5, 3, 4 and 5mm
- 1 Spare bulb, 225 V, 15 Watt
- 1 Grease gun with grease
- 1 T-piece, for vacuum connection
- 1 Plastic tube, 1/2" x 1500mm with hose clamps
- 1 Plastic tube, 3/4" x 1500mm with hose clamps
- 1 Armored plastic tube 3/8" x 1500mm with hose clamp for vacuum pump
- 1 Dial gauge with holder for lining up diamond cup wheel
- 1 Holder for specimens of irregular shapes (max. 75x75mm)
- 1 Movable stop
- 1 Vacuum chuck for resectioning
- 1 Splash guard for the cutting module
- 6 Polyethylene washers for the diamond cup wheel
- 1 I ADDUN additive
- 3 Screws M6 x 20 for mounting of diamond cup wheel

1.2 Mounting and Assembly

1.2.1 Discoplan-TS

- Remove the "Tectyl" protecting film with kerosene or petrol.
- Place Discoplan-TS on top of a sturdy table. Make sure that the machine is in level to avoid warping.

1.2.2 Recirculating Cooling Unit

- Fill the recirculating tank with water, minimum 180mm high.
- Add cooling liquid additive. It is important to follow the instructions for use on the additive bottle. Shake the bottle well before the additive is added to the water.
- Place the lid with the pump on the recirculating tank.
- Mount the hoses of the recirculating cooling unit on the pipe branches (fig. 4.7 and 4.8).
- Lead the hoses to the branches of the recirculating cooling unit. Avoid sharp bends on the hoses.

1.3 Electrical Connection

The cable is connected to Discoplan-TS in the connection box (fig. 4.9). See diagram in the lid of the connection box.

Note: The single phase version of Discoplan-TS (TISEN) is connected with a Recirculating Cooling Unit (TRECI) for three phases. The required capacitor is built into Discoplan-TS.

The mains cable (fig. 4.16) is connected to a plug switch.

IMPORTANT

Before Discoplan-TS is connected, check that the mains voltage corresponds to the voltage indicated on the type plate on the rear of the machine.

1.4 Vacuum Connections

Connect the vacuum chuck for resectioning to the branch (fig. 4.20). Discoplan-TS obtains the vacuum effect through a pump with armored plastic tubes connected to the branches (fig. 4.18 and 4.19).

Remember to tighten the hose clamps!

2. Operation

Table of Contents	Page
2.1 Introduction	8
2.2 Controls	9
2.3 The TS-Method for Preparing Thin Sections.	10
2.3.1 Necessary Struers Equipment	
2.3.2 The TS-Method Steps	
2.4 Mounting of Wheels.	11
2.4.1 Mounting the Diamond Cut-off Wheel	
2.4.2 Mounting the Diamond Cup Wheel	
2.5 General Operation of Discoplan-TS.	12
2.5.1 Start of Motor	
2.5.2 Cooling Water	
2.5.3 Cutting a Specimen	
2.5.4 Grinding a Glass Slide	
2.5.5 Resectioning	
2.5.6 Grinding a Mounted Specimen	

2.1 Introduction

Discoplan-TS is a mineralogical cutting and grinding machine for fast and accurate production of thin sections. With Discoplan-TS it is possible to carry out a very precise grinding of glass slides and thin sections.

Precision is obtained through an absolutely stable construction. All vital parts of Discoplan-TS are made of heavy vibration-dampening cast iron. The diamond wheels are driven by a powerful motor, especially designed for vibration-free and smooth running.

2.1.4 Cutting Module

A diamond cut-off wheel is mounted in the left module of Discoplan-TS. The cutting module has a cutting table with guide ways and guide bar, used during the cutting process. Standard equipment includes a holder for cutting rocks and a movable stop for cutting standard specimens (8 x 20 x 30mm). A separate vacuum chuck is used for resectioning.

2.1.2 Grinding Module

A diamond cup wheel for precision grinding is mounted in the right module of Discoplan-TS. The grinding module has 3 vacuum chucks for grinding of glass slides and thin sections. The chucks are mounted on a feed shaft controlled by a micrometer screw.

2.1.3 Recirculating cooling unit

Discoplan-TS is delivered with an effective recirculating cooling unit. The unit consists of a 35 l tank and a lid with a sedimentation system. The strong 110 watt pump has adequate cooling capacity for both cutting and grinding. Discoplan-TS has individual adjusting knobs for the cutting and grinding module cooling (fig. 1.6 and 1.7).

2.2 Controls

2.2.1 Adjusting Knobs

Fig. 1.6

Cooling water for cutting module: adjusting knob for the cutting module.

Fig. 1.12

Motor: the motor switch turns the central motor on and off.

Fig. 1.7

Cooling water for grinding module: adjusting knob for the grinding module.


2.2.2 Other Controls

Fig. 1.3

Micrometer screw: control the grinding depth of the diamond cup wheel in the grinding module.

Fig. 1.9

Vacuum button: release the thin section in the vacuum chuck of the cutting module.

2.3 The TS-Method for Preparing Thin Sections

2.3.1 Necessary Struers Equipment

The TS-Method (Thin Section Method) developed by Struers involves the following Struers equipment:

- A Discoplan-TS mineralogical cutting and grinding machine (the machine described in this Instruction Manual)
- Glass slides and mounting resin for the specimens¹⁾
- A Planopol-V²⁾ lapping/polishing machine with a PdM-Force specimen mover.
- Lapping tools comprising cast iron disc and lapping powder³⁾
- Thin section holders with boron carbide sticks (BORTY)⁴⁾ and thin section holders for polishing (TYNDS)⁵⁾
- Epovac⁶⁾ vacuum impregnation apparatus
- Drybox-2⁷⁾ specimen drier
- Glass slides⁸⁾

¹⁾ To be purchased separately - see Section 3, Equipment, Accessories and Consumables

2.3.2 The TS-Method Steps

- ① Sampling on Discoplan-TS¹⁾
- ② Cutting of specimens 8x20x30mm on Discoplan-TS.
- ③ Automatic lapping of one side of specimen on Planopol-V with PdM-Force.
- ④ Glass slide.
- ⑤ Grinding of glass slide in Discoplan-TS to a known thickness x, (e.g. x = 1.164mm).
- ⑥ Mounting lapped side of specimen onto the ground glass slide in Struers vacuum impregnation apparatus Epovac.
- ⑦ Cutting off surplus material of specimen on Discoplan-TS (e.g. t' = 0.5 - 1mm).
- ⑧ Grinding of thin section in Discoplan-TS to the thickness t' including glue (e.g. t' = 0.080mm).
- ⑨ Automatic lapping of thin section on Planopol-V with PdM-Force and Thin Section Holder with boron carbide sticks. Lapping to thickness t'' including glue (e.g. 1.204 = x + t'' + z, where z is 10µm, which is the amount removed by the abrasive particles below the level of the boron carbide sticks with 1000 grit, z is normally approx. t'' = 1.204 - 1.164 - 0.010 = 0.030mm = 30µm).
- ⑩ Automatic polishing of thin section in thin section holder on Planopol-V with PdM-Force or mounting of cover glass²⁾. During polishing approx. 10µm are removed (e.g. t''' = t'' = 20µm).

The finished specimen (11) is now 20µm, provided the stated figures are used (the 5µm Epofix glue layer is determined experimentally).

¹⁾ Vacuum impregnation in Struers Epovac may be necessary. Use Struers Epofix for mounting.

²⁾ Polishing may be omitted and the specimen covered with a cover glass instead. The required end thickness t''' (specimen + cement) in point 9 must then be approx. 30µm.

2.4 Mounting the Wheels

2.4.1 Mounting the Diamond Cut-off Wheel

- Loosen the thumb screws (fig. 5.10) and remove the cutting table (fig. 1.2).
- Remove nut, washer and flange and mount the cut-off wheel, flange and washer.
- Tighten the nut with the keys supplied with Discoplan-TS (fig. 5.9 and 5.12) and mount the cutting table.

2.4.2 Mounting the Diamond Cup Wheel

- Tilt the 3 vacuum chucks (fig. 1.4) outward with the handle and dismount the splash guard (fig. 1.8).
- Remove the screws and washers from the flange (fig. 6.1) with the 5mm Allen key supplied with Discoplan-TS.
- Introduce the screws through the holes in the diamond grinding wheel, and put 3 polyethylene washers¹ round the screws from the rear side.
- Place the diamond cup wheel on the flange (fig. 6.1) and tighten it.
- Clamp the holder (fig. 7.3) for the dial gauge (fig. 7.1) on the protection guard (fig. 7.2).
- Place the dial gauge in the holder (the placing of the dial gauge can be modified vertically).
- Adjust the diamond cup wheel by tightening the screws and by compressing the 3 polyethylene washers. This has to be done gradually.

IMPORTANT

- The adjustment should be carried out meticulously and within 1/100mm. Struers recommends that you check the adjustment regularly.
- When the diamond cup wheel has been adjusted, the vacuum chucks should be plane ground.

¹ The washers must be of equal thickness which is ensured by measuring and sorting out the washers.

2.5 General Operation of Discoplan-TS

2.5.1 Start of Motor

- Turn the MOTOR knob to ON. The motor driving both the cutting and grinding module starts.
- Turn both cooling water knobs on and let Discoplan-TS run for about half an hour. This ensures that all parts of the machine have obtained thermal equilibrium.

2.5.2 Cooling Water

- Adjust the position of the left (cutting module) and right (grinding module) cooling water tubes manually to fit the size and position of the specimen.
- Adjust the water pressure on the cooling knobs to the left (cutting module) and the right (the grinding module). Make sure that the water tap fig. 6.2 points to the diamond cup-wheel.

2.5.3 Cutting the Specimen

- Loosen the guide bar thumb screws and make room for the specimen. Place the specimen pointing at the diamond cut-off wheel. Adjust the position of the specimen by moving the guide bar and the holder.
- Fasten the guide bar and the holder with the thumb screws.
- Cut the specimen by manually pressing it against the cut-off wheel.
Always wear safety goggles when working with Discoplan-TS.

Cutting in General

Large Rocks

Large rocks (max. 75x75mm) are either led by hand or in the holder (fig. 8.1). The holder is guided by the adjustable guide bar (fig. 8.2).

Standard specimens (8x20x30mm)

After having cut out an adequate specimen of the rock, standard specimens can be cut using the movable stop (fig. 9.1):

- Hold the specimen against the inner corner of the movable stop.
- Cut the specimen by moving the stop in the feed track. The track to the far right controls the length, the middle track controls the width and the left track controls the depth of the specimen.
- The specimens are now ready for lapping (see the TS-Method description).

2.5.4 Grinding the Glass Slides

- Make sure that the vacuum chuck surface is plane. The vacuum chucks should be dressed regularly at each 20 glass slide grinding session (see Section Maintenance, Dressing of Vacuum chucks).
- Check the planeness of the glass slide with a straight rail. Unplane glass slides should be scrapped.
- Lower the handle to open the splash guard. Mount the glass slides on the vacuum chucks. Raise the handle again to close the splash guard.
- Move the handle lightly (use 2-3 fingers and do not press hard!) while adjusting the micrometer screw in small steps until you can hear the characteristic sound of the glass slides against the diamond cup wheel.
- Grind at a steady and even feed in order to permit the diamond cup wheel to cut the glass slide.
- Do not move the micrometer screw more than $\leq 0.01\text{mm}$ because of deflections in the machine.
- Grind the slides 10-20 times to ensure complete planeness.
- Adjust the micrometer screw again and grind another 10 times.
- Repeat the adjusting/grinding until the desired glass thickness has been reached.
- After grinding the glass slide to the desired thickness, set the micrometer scale to zero. A glass slide sample has been enclosed with Discoplan-TS for your convenience. The glass slide sample enables you to double-check the correct thickness of your glass slides.
- Remove the glass slides by pressing the vacuum button
- Cement the lapped standard specimen to the ground face of the glass slide (see the TS-Method). Allow the cement to harden, and carry out the resectioning on the cutting module of Discoplan-TS.

2.5.5 Resectioning

- Mount the glass slide with the specimen in the vacuum chuck. Place and adjust the guide-bar (fig. 10.2).
- Cut the vacuum chuck with the specimen against the rotating diamond cut-off wheel.
- After the cutting process remove the glass slide with the specimen by pushing the vacuum button (fig. 10.3). The specimen is now ready for grinding in the grinding module (See the TS-Method).

2.5.6 Grinding a Thin Section

- As the micrometer scale has been set to zero after the grinding of the glass slide it is possible to read the final thickness of cement + specimen directly on the scale during the grinding of the thin section (normally 80µm).
- Lower the handle to open the splash guard and mount the thin section on the vacuum chucks.
- Turn the micrometer screw 5-10 turns away from the diamond cup wheel. Count the number of turns and write down the total. Raise the handle again to close the splash guard.
- Move the handle lightly (use 2-3 fingers and do not press hard!) while adjusting the micrometer screw in small steps until you can hear the characteristic sound of the thin section against the diamond cup wheel.
- Grind at a steady and even feed in order to permit the diamond cup wheel to cut the specimen.
- As the thickness of the thin section decreases you decrease the micrometer screw steps accordingly until you are close to the final thickness.
- When you get close to the final thickness of the specimen, you should grind the thin section 10-20 times per step to ensure complete planeness. Do not move the micrometer screw more than ≤ 0.01 mm because of deflections in the machine.
- When you have removed the material to the final thickness of the thin section the micrometer screw should be back in the zero position from the grinding of the glass slide plus 80µm.
- Remove the thin section by pressing the vacuum push button.
- Proceed with lapping and polishing on Planopol-V according to the TS-Method.

3. Equipment, Accessories and Consumables

3.1 Accessories and Consumables for Discoplan-TS

The cutting module takes diamond cut-off wheels in the range of $\varnothing 150\text{mm}$ to $\varnothing 230\text{mm}$, with a thickness from 0.5 - 1.5mm. Struers, however, recommends a diamond cut-off wheel $\varnothing 200\text{mm} \times 1\text{mm}$ in order to obtain the required flexibility. Center hole 22mm.

The $\varnothing 200 \times 1 \times 22\text{mm}$ is a quality cut-off wheel with the diamond grains bonded in metal along the wheel periphery.

For the grinding module we recommend a metal-bonded diamond cup wheel (max. $\varnothing 180\text{mm}$). Struers diamond cup wheel (01TYN), $\varnothing 150\text{mm}$, metal bonded, $70\mu\text{m}$, specially made for Discoplan-TS, has good free-cutting qualities.

A $35\mu\text{m}$ bakelite bond diamond cup wheel (06TYN) is also available.

Specification	Code
Accessories	
<i>Laboratory Furniture</i> Table unit. With a rolling pallet for recirculating cooling unit	LABDI
<i>Vacuum Pump</i> used with vacuum chucks on Discoplan-TS. With vacuum: 780 mbar and output: 5.8 l/min.	TISPU
Consumables	
<i>Diamond Cut-off Wheel</i> for cutting of minerals and ceramics (> HV 800). Metal bond. $\varnothing 200\text{mm} \times 1.0\text{mm} \times 22\text{mm}$	20TRE
<i>Diamond Cut-off wheel</i> for cutting of sintered carbides and ceramics (> HV 800). Bakelite bond. $\varnothing 200\text{mm} \times 1.1\text{mm} \times 22\text{mm}$	23TRE
<i>Diamond Cup Wheel, $70\mu\text{m}$</i> Metal bonded, for grinding of hard, brittle materials on Discoplan-TS. $\varnothing 150\text{mm}$	01TYN
<i>Diamond Cup Wheel, $35\mu\text{m}$</i> Resin bonded, for grinding of hard, ductile materials on Discoplan-TS. $\varnothing 150\text{mm}$	06TYN

3.2 Equipment and Accessories for the TS-Method

The following equipment is necessary to take full advantage of the TS-method.

Specification	Code
Vacuum Impregnation	
<i>Epovac</i> vacuum impregnation apparatus with air ejector pump. Complete with vacuum chamber, pressing rod for gluing, and consumables for impregnation. Compressed air supply 4-8 bar required	EPOMA
<i>Epovac</i> vacuum impregnation apparatus. Complete with vacuum chamber, pressing rod for gluing, and consumables for impregnation. External vacuum pump required	EPOCA
Specimen drier	
<i>Drybox-2</i> combined specimen drier and curing oven with adjustable temperature (approx. 40-90°C) and low noise fan. With basket for specimens and slot for specimen holders	DESIC
Grinding/polishing machine	
<i>Planopol-V</i> universal grinding/polishing machine with variable speed (25-800 rpm). For disc 300mm dia. With water inlet and outlet, splash ring and water tap.	PLAVE
Specimen Mover	
<i>PdM-Force</i> Automatic specimen mover for mineralogical specimens. For mounting on Planopol-V. Lubricator included.	FORCE
<i>Specimen holder plates</i> for ø300mm polishing disc with holes for 7 petrographic specimens 20 x 30mm for 5 specimens ø60mm or 5 Thin Section Holders.	FORON FORIX
<i>Weights</i> set with pressure rod 1 set, 6 x 200 g	FORTO

Specification	Code
Thin Section Holders	
<i>Thin Section Holder with Boron Carbide Sticks</i> for lapping of thin sections, ø60mm. Available for standard glass slides 27 x 46mm, 28 x 48mm and 30 x 45mm. Including 4 foils of 25µm and 3 foils of 8µm	BORTY
<i>Thin Section Holder</i> for polishing and lapping of thin sections, ø60mm. Available for standard glass slides 27 x 46mm, 28 x 48mm and 30 x 45mm	TYNDS

3.3 Consumables for the TS-Method

Specification	Code
Epovac	
<i>Epofix Kit</i> consisting of 1 l resin, 130 ml hardener and consumables	EPOFI
<i>Caldofix Kit</i> consisting of 1 l resin, 400 ml hardener and consumables	CALDO
Drybox-2	
<i>Lakeside 70 C Cement</i> natural thermoplastic resin. Fluid at 140°C. Package of 12 bars	CEMEN
<i>Canada Balsam</i> natural thermoplastic resin, pure. 100 ml	BALHU
<i>Eukitt</i> synthetic thermoplastic resin for mounting of cover glasses on thin sections. 250 ml.	EUKIT
Planopol-V	
<i>Cast Iron Lapping Discs</i> concentrically grooved disc, for lapping of mineralogical specimens. ø300mm ø200mm	GRICO GRIAL
<i>Silicon Carbide Powder</i> for lapping of mineralogical specimens. To be used on cast iron lapping disc. Grain sizes: 120, 220, 320, 400, 600, 800, 1000, 1200. 500 g (grain 120-800) and 400 g (grain 100 and 1200)	POWDE

Specification	Code
<i>Aluminium Oxide Powder</i> for lapping of mineralogical specimens. To be used on cast iron lapping disc. Grain sizes: 120, 220, 320, 400, 500 g	POWAL
<i>Petrodisc</i> polishing disc for plane polishing of mineralogical and ceramic specimens ø290mm ø175mm	PETRO PETOW
<i>Dressing Stick</i> Aluminium oxide stick for trueing and opening of diamond and CBN cut-off wheels, diamond grinding discs, Diamond Pad and Petrodisc-M	OPRET
<i>DP-Spray, P</i> high performance diamond product containing exclusively polycrystalline diamonds, in spray can grain size 6µm grain size 3µm grain size 1µm <i>DP-Spray, M</i> High performance diamond product containing monocrystalline diamonds, in spray can grain size 6µm grain size 3µm grain size 1µm	SPRIX SPRET SPRON SMAIX SMAET SMAON
Glass Slides	
<i>Standard Slides</i> for thin sections, with ground edges, 27 x 46mm, 28 x 48mm and 30 x 45mm. 100 pcs.	OBJEK
<i>Microcover Glass</i> to be used for covering thin sections. 24 x 32mm. Approx. 100 pcs.	COVER

4. Maintenance

4.1 Daily (8 Hours of Use)

4.1.1 Dressing the Diamond Wheels

Dress the diamond wheels with a dressing stick (aluminium oxide) in order to ensure that the wheels always have the same good cutting properties.

4.1.2 Plane Grinding of the Vacuum Chucks (Grinding Module)

The ceramic vacuum chucks must be dressed after adjustment of the diamond cup wheel or before grinding of plan parallel glass slides.

- Draw a cross on the surface of the vacuum chucks with a pen.
- Carefully grind the vacuum chucks with the diamond cup wheel without stressing neither the handle, the micrometer shaft nor the machine itself.
- Do not remove more than 10-20µm at a time as a maximum.
- After the last turn of the micrometer screw, you should grind at least 10 times.
- The truing level is correct when the whole cross drawn with the pen has disappeared and a parallelism has been obtained.

4.2 Weekly (40 Hours of Use)

Discoplan-TS and the recirculating cooling unit should be cleaned regularly.

Discoplan-TS must be kept free from rust and it is important to follow the instructions on the bottle with cooling liquid additive.

Grease the feed shaft with the grease gun.

4.3 Monthly (160 Hours of Use)

The bearings (fig. 6.5) can be readjusted, if needed, with the Allen key supplied with Discoplan-TS, width 5mm. The micrometer shaft must not run too tight in the bearings.

4.4 Replacements

Bulb

The bulb (15 W, fig. 11.14) can be replaced after the cover plate (fig. 5.3) has been removed. This is done by loosening the screw (fig. 4.13) and drawing the plate backwards and upwards.

Reflector

The same applies to the reflector (fig. 11.16), which can be adjusted or removed for cleaning after the screw (fig. 11.15) has been loosened (use the 3mm Allen key supplied with Discoplan-TS).

Micrometer Screw

Adjusting axial play in the micrometer screw is done by tightening the end nut after removing the cover (fig. 6.6).

5. Technical Data

5.1 Discoplan-TS

<i>Voltage</i>	3 x 220V, 50-60Hz 3 x 380V, 50Hz 3 x 415V, 50Hz 3 x 440V, 60Hz 1 x 220V, 50Hz
<i>Power</i>	736W (3-phase) 442W (1-phase)
<i>Rotational Speed</i>	1400 rpm at 50 Hz 1700 rpm at 60 Hz
<i>Diamond cut-off wheel</i>	ø150mm - ø230mm Thickness: 0.5-1.5mm (0.02"-0.06")
<i>Diamond cup wheel</i>	Max. ø180mm
<i>Micrometer Screw</i>	Precision scale including vernier: 1µm
<i>Vacuum chucks</i>	3 chucks 30 x 50mm
<i>Weight</i>	65 kg
<i>Dimensions</i>	Height: 320mm Width: 700mm Depth: 370mm

5.2 Recirculation Cooling Unit

<i>Voltage</i>	3-phase
<i>Power</i>	140 Watt
<i>Dimensions</i>	Height: 410mm Width: 500mm Depth: 400mm
<i>Tank Volume</i>	35 l
<i>Weight</i>	7 kg

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Drawing

Fig. 1
Fig. 2
Fig. 3
Fig. 4
Fig. 5
Fig. 6
Fig. 7
Fig. 8
Fig. 9
Fig. 10
Fig. 11

1.1	Discoplan-TS
1.2	Front Panel
1.3	Recirculation Cooling Unit
1.4	Back
1.5	Cutting Module
1.6	Grinding Module
1.7	Micrometer Screw
1.8	Cutting Module, Details
1.9	Cutting Module, Holder
1.10	Vacuum Holder, left
1.11	Inside View
1.12	Various

The drawings are not to scale. Some of the drawings may contain position numbers not used in connection with this manual.


Drawing

Fig. 1

Pos. Description

Spare Part No.

1.1 Front Plate

Drawing

Fig. 2

Pos.	Description	Spare Part No.
1.2	Front Panel	
	Arm, complete with valves and ceramic plates	372MP001
	Front splash guard	372MP035
	Set of vacuum hoses (external)	372MP087

Drawing

Pos. Description

Spare Part No.

1.3 Recirculation Cooling Unit

See accompanying Manual (TRECI)

Drawing


Fig. 4

Pos. Description

Spare Part No.

1.4 Back


17 Handle372MP084

Drawing

Fig. 5

Pos.	Description	Spare Part No.
1.5 Cutting Module		
9	Hook-spanner for mounting of diamond cut-off wheel	372MP040
10	Thumb screw	372MP088
10	Cutting table	13729038
12	Spanner for mounting of diamond cut-off wheel	372MP041
13	Outside flange	372MP037
	Nut and disk	372MP038
13	Inside flange	372MP039
14	Motor (state voltage and cycles)	372MP003
	Ring for oil seal, 2 pcs. (placed in each end of the motor spindle)	372MP093
	Ring for oil seal, 2. pcs. (placed in each end of the motor spindle)	372MP93
15	Axle	372MP036



Drawing

Fig. 6

Pos. Description
1.6 Grinding Module

Spare Part No.

Screws for diamond cup wheel 372MP091
 Polyethylene washers (6 pcs.) 372MP092


Drawing

Fig. 6A

Pos.	Description	Spare Part No.
Grinding Module		
1	Turntable scale.372MP042
2	Spindle guide.372MP043
3	Hexagon-headed screw MC4x20372MP044
4	Bearing frame, right	13729034
85	Bearing frame, left	13729043
5	Retaining ring372MP046
6	Pressure spring372MP047
8	O-ring ø3.3x2.4.372MP049
9	Value plate372MP050
	3 Hexagon-headed screw MSF 3x5	
	4 Hexagon-headed screw MC 6x20	
10	Arm for vacuum holder.372MP051
	2 Hexagon headed screws MC 6x16	
11	Spring retainer372MP052
12	Protecting hose.372MP053
13	Spring retainer(=372MP052)372MP054
14	Sealing ring.372MP055
15	Bearing372MP056
16	Cover for bearing housing372MP057
17	Bearing ring372MP058
18	Axle nut.372MP059
19	Bearing housing372MP060
20	Spindle tube372MP061
21	Hexagon-headed screws (3 pcs.)372MP062
22	Spindle with bronze bearing and safety nut.372MP063
23	Bronze bearing.372MP064
24	Intermediate disk, bronze372MP065
25	Clamping screw372MP066
26	Thumb screw.372MP067
27	Sealing ring (=372MP055).372MP068
28	Retaining ring (=372MP046).372MP069
29	Stop disk372MP070
30	O-ring ø3,53x29,75372MP071
31	Spring (2 pcs.)372MP072
32	Cylindrical pin ø8x18.372MP073
33	Cylindrical pin ø6x30.372MP074
34	Lever arm with handle372MP075
37	Hexagon-headed screw (2 pcs.)372MP076
	Grease gun with grease372MP090
38	Cylindrical pin ø2x14.	2ZS01214

Drawing 6A

Drawing


Fig. 6B

Pos. Description

Spare Part No.

Grinding Module

1	Foot-lever.....	12480410
2	Ceramic vacuum plate with 4 head-screws M4x16, stainless 4 Nylon discs ø7xø4x1 High vacuum grease.....	372MP078 (TISPA)
3	L-Fitting 1500 2-2.....	372MP079
4	Valve spindle, long.....	372MP080
5	Valve spindle, short.....	13728011
6	Hose fitting 11752-1.....	13728049
	Vacuum holder complete, including arm.....	372MP001

Drawing


Fig. 7

Pos. Description

Spare Part No.

1.7 Micrometer Screw


1	Dial gauge398MP007
3	Holder for dial gauge, complete372MP026

Drawing

Fig. 8


Pos.	Description	Spare Part No.
1.8 Cutting Module, details		
1	Holder for large rocks372MP033
2	Guide plate372MP032
7	Safety guard, left side372MP029
9	Thumb screw372MP028

Drawing

Fig. 9

Pos.	Description	Spare Part No.
1.9	Cutting Module, Holder	
1	Movable stop.372MP034

Drawing


Fig. 10

Pos. Description

Spare Part No.

1.10 Vacuum Holder, Left

1	Vacuum Holder, complete (TISOM)	03726902
3	Valve spindle.....	13728011
	O-rings for valve spindle	372MP049
4	Vacuum chuck	13720051
5	Guard fixture	13140281


Drawing

Fig. 11

Pos.	Description	Spare Part No.
1.11	Inside view	
1	Splash guard, top.....	372MP024
2	Splash guard, outer.....	372MP086
4	Contact, state voltage and cycles.....	372MP004
4	Thermal switch, state voltage and cycles.....	372MP089
5	Switch	372MP005
6	Water cock, Cutting Module.....	372MP006
7	Water cock, Grinding Module.....	372MP006
8	Inlet tube, complete.....	372MP012
9	Outlet tube, complete.....	372MP013
	Inlet hose	372MP014
	Outlet hose.....	372MP015
	Set of internal hoses	372MP016
10	Tube for cooling water, left side	372MP027
11	Water nozzle, right side.....	372MP019
13	Circuit board, complete (state motor voltage and cycles) ..	372MP010
14	Filament lamp.....	372MP009
16	Light conductor.....	372MP011
17	Arm for microswitch.....	372MP018
18	Microswitch with arm.....	372MP008
19	Adjusting screw, complete.....	372MP020
20	Vacuum tube, complete.....	372MP030
22	Vacuum tube, complete.....	372MP021
	Mounting flange for diamond cup wheel, complete.....	372MP022
	Spring for splash guard	372MP023
23	Vessel for cooling water	372MP031

Inside View

Drawing

Pos. Description

Spare Part No.

1.12 Various


Vacuum pump TISPU (230V / 50Hz)	03726433
Vacuum pump TISPU (115V / 60Hz)	03726421
Filter for vacuum pump	13722901
Vacuum hose. PVC tube, clear, 3/16" (10 M)	RNU19306

2. Diagrams

Table of Contents

Drawing

2.1	Electrical diagram, 3-phase	3720095
2.2	Electrical diagram, 1x220V/50 and 60Hz.....	12480527
2.3	Electrical diagram, 1x115V/60Hz.....	12480528

Matr.:

Erstatning af 225V-15W
509. 29.9.86 770/14

AB. Blum + S. R. 273 + afkorting af aksel bitt
210696 877/77

STRUERS

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Målestok

Tegn.


MG/BS

Dato

24.3.71

Diagram, Discoplan
1 x 220V/50Hz og 60Hz

12480527A

Matr.:

STRUERS

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Målestok

Tegn.

HSB/RP

Dato

12.10.78

Diagram, Discoplan
1x 115 V/60Hz

12480528A

Erstatte 3V4-M-5/0
29.9.86 7803 / 2/4

A: brun, sort 3+ sort 2 til f. 2/6 9 BHT/4V