

# **ACQUITY UPLC H-Class and H-Class Bio**

## **System Specifications**

Revision B

**Waters**  
THE SCIENCE OF WHAT'S POSSIBLE.™

Copyright © Waters Corporation 2010  
All rights reserved

## Copyright notice

---

© 2010 WATERS CORPORATION. PRINTED IN THE UNITED STATES OF AMERICA AND IN IRELAND. ALL RIGHTS RESERVED. THIS DOCUMENT OR PARTS THEREOF MAY NOT BE REPRODUCED IN ANY FORM WITHOUT THE WRITTEN PERMISSION OF THE PUBLISHER.

The information in this document is subject to change without notice and should not be construed as a commitment by Waters Corporation. Waters Corporation assumes no responsibility for any errors that may appear in this document. This document is believed to be complete and accurate at the time of publication. In no event shall Waters Corporation be liable for incidental or consequential damages in connection with, or arising from, its use.

# Table of Contents

---

|  |    |
|--|----|
| Copyright notice .....                 | ii |
| Introduction .....                     | 1  |
| System features .....                  | 2  |
| Instrument control .....               | 4  |
| Environmental specifications .....     | 4  |
| Electrical specifications .....        | 5  |
| Performance specifications .....       | 6  |
| Quaternary solvent manager .....       | 6  |
| Sample manager - FTN .....             | 9  |
| Column heater .....                    | 16 |
| Column manager.....                    | 18 |
| Wetted materials of construction ..... | 21 |
| Quaternary solvent manager .....       | 21 |
| Sample manager - FTN .....             | 21 |
| Column heater .....                    | 22 |
| Column manager.....                    | 22 |



# Introduction

---

The system specifications outlined in this document depend on the conditions in individual laboratories. Refer to the *ACQUITY UPLC H-Class Site Preparation Guide*, or contact Waters Technical Service for more information on specifications.

## Notes:

- If your system includes an ELS detector, see the *ACQUITY UPLC Evaporative Light Scattering Detector Getting Started Guide* for physical specifications.
- If your system includes a PDA detector, see the *ACQUITY UPLC Photodiode Array Detector Getting Started Guide* for physical specifications.
- If your system includes a PDA Extended  $\lambda$  detector, see the *ACQUITY UPLC Photodiode Array Extended  $\lambda$  Detector Getting Started Guide* for physical specifications.
- If your system includes an FLR detector, see the *ACQUITY UPLC Fluorescence Detector Getting Started Guide* for physical specifications.
- If your system includes a mass spectrometer, see the documentation included with it for specifications.

## System features

---

The following table lists the ACQUITY UPLC® H-Class and H-Class Bio system features.

### System features

| Item                       | Specification   |
|----------------------------|---|
| Dwell (delay) volume       | <400 µL with 100-µL mixer   |
| Integrated leak management | Drip trays direct all leaks to the front of the instrument and then into waste line.  |
| Leak detection             | Leak sensors, installed in drip trays.  |
| Quantum synchronization    | Injection synchronization between pump and sample manager. Enhances retention time reproducibility.   |
| Settable flow rate range   | 0.010 to 2.000 mL/min, in 0.001-mL increments   |
| Maximum operating pressure | 103,421 kPa (1034 bar, 15,000 psi) to 1 mL/min;<br>62,053 kPa (621 bar, 9000 psi) to 2 mL/min   |
| pH range                   | 2 to 12   |
| Cycle time                 | <p>&lt;30 s</p> <p>System cycle time (or overhead) is equal to the chromatographic run time subtracted from the injection-to-injection time.</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• System: ACQUITY UPLC H-Class quaternary solvent manager (QSM), ACQUITY UPLC H-Class sample manager with flow through needle (SM-FTN), ACQUITY UPLC H-Class column heater with active pre-heater (CH-A), and ACQUITY UPLC TUV detector</li> <li>• Isocratic chromatography</li> <li>• Flow rate: ≥0.4 mL/min</li> <li>• Injection volume: 1 µL</li> <li>• SM-FTN parameters: Default aspiration speeds and wash times</li> </ul> |

## System features (Continued)

| Item                     | Specification  |
|--------------------------|--|
| Cycle time (continued)   | <ul style="list-style-type: none"> <li>• Load ahead mode: enabled</li> <li>• Loop offline: 0.2 min</li> <li>• Run time: 2.0 min</li> </ul>   |
| Gradient mixers          | <p>H-Class QSM:</p> <ul style="list-style-type: none"> <li>• Standard: stainless steel, 100-μL mixer/filter</li> <li>• Optional: stainless steel, 250-μL mixer/filter</li> </ul> <p>H-Class bioQSM:</p> <ul style="list-style-type: none"> <li>• Standard: titanium, 100-μL mixer/filter</li> <li>• Optional: titanium, 250-μL mixer/filter</li> </ul>   |
| Plunger wash feature     | Wash pump plungers using seal-wash solvent, can be primed manually or run automatically.   |
| No-flow shutdown feature | Automatically runs the wash plungers function after a user-specified period of idle time.  |
| Column tracking          | eCord™ Technology column information management tracks and archives column usage history   |
| Unattended operation     | Leak sensors, full 96-hour diagnostic data displayed through ACQUITY UPLC console software   |
| Auto•Blend Plus™         | Software with which users can blend liquid mobile phases (acid, base, aqueous, and salt solution) of known composition and concentration. Users select individual steps with a specified pH and/or salt concentration at a particular time and enter one of 11 predefined functions to transition between these times and conditions. The Auto•Blend Plus technology then automatically calculates the percent acid, percent base, percent salt, and percent solvent required to deliver the specified pH and salt conditions at a given time. |

## Instrument control

---

The following table lists the mechanisms used to control ACQUITY UPLC H-Class and H-Class Bio system instruments.

### Instrument control

| Item                    | Specification  |
|-------------------------|--|
| External control        | Empower™ software, MassLynx™ software, or standalone, through ACQUITY UPLC Console software  |
| External communications | Ethernet interfacing via RJ45 connection to host PC  |
| Event inputs/outputs    | Rear panel contact closure and/or TTL inputs/outputs   |
| Connections INSIGHT®    | Provides real-time monitoring and automatic notification of instrument performance and diagnostic information, allowing for quicker problem resolution |
| Local control           | ACQUITY UPLC Local Console Controller (LCC)  |

## Environmental specifications

---

The following table lists the environmental specifications for the ACQUITY UPLC H-Class and H-Class Bio instruments.

### Environmental specifications


| Attribute                        | Specification                |
|----------------------------------|------------------------------|
| Acoustic noise                   | <65 dBA, system              |
| Operating temperature            | 4 to 40 °C (39.2 to 104 °F)  |
| Operating humidity               | 20 to 80%, noncondensing     |
| Shipping and storage temperature | –30 to 60 °C (–22 to 140 °F) |
| Shipping and storage humidity    | 20 to 80%, noncondensing     |



## Electrical specifications

The following table lists the electrical specifications for the ACQUITY UPLC H-Class and H-Class Bio instruments.

### Electrical specifications

| Attribute  | Specification   |
|--|---|
| Protection class <sup>1</sup>  | Class I   |
| Overvoltage category <sup>2</sup>  | II  |
| Pollution degree <sup>3</sup>  | 2   |
| Moisture protection <sup>4</sup>   | Normal (IPXO)   |
|  Line voltages, nominal | Grounded AC   |
| Voltage range  | 90 to 264 Vac   |
| Frequency  | 47 to 63 Hz   |
| Maximum power draw   | QSM: 360 VA<br>SM-FTN: 400 VA<br>Column manager with active pre-heater (CM-A): 400 VA |

- Protection Class I** – The insulating scheme used in the instrument to protect from electrical shock. Class I identifies a single level of insulation between live parts (wires) and exposed conductive parts (metal panels), in which the exposed conductive parts are connected to a grounding system. In turn, this grounding system is connected to the third pin (ground pin) on the electrical power cord plug.
- Overvoltage Category II** – Pertains to instruments that receive their electrical power from a local level such as an electrical wall outlet.
- Pollution Degree 2** – A measure of pollution on electrical circuits that can produce a reduction of dielectric strength or surface resistivity. Degree 2 refers only to normally nonconductive pollution. Occasionally, however, expect a temporary conductivity caused by condensation.
- Moisture Protection** – Normal (IPXO) – IPXO means that no Ingress Protection against any type of dripping or sprayed water exists. The “X” is a placeholder that identifies protection against dust, if applicable.

## Performance specifications

---

### Quaternary solvent manager

The following table lists the performance specifications for the ACQUITY UPLC H-Class QSM and H-Class bioQSM.

#### QSM performance specifications

| Item                 | Specification  |
|----------------------|--|
| Number of solvents   | One to four (A,B, C, and D), in any combination. Optional 6-position solvent selection valve enables solvent selections D <sub>1</sub> through D <sub>6</sub> on line D, in addition to A, B, and C (a total of nine solvents to select from).   |
| Solvent conditioning | Integrated vacuum degassing, four chambers. One additional for the SM-FTN.   |
| Gradient formation   | Low-pressure mixing, quaternary gradient   |
| Gradient profiles    | 11 gradient curves, including linear, step (2), concave (4), and convex (4)  |
| Primary check valve  | Intelligent Intake Valve ( <i>i<sup>2</sup>Valve</i> )   |
| Flow accuracy        | ±1.0% of set flow at 0.500 to 2.000 mL/min using 100% solvent A. Backpressure 4137 to 6895 kPa (41 to 69 bar, 600 to 1000 psi), with water.  |
| Flow precision       | 0.075% RSD or ±0.02 min SD, whichever is greater, based on six replicates.<br>Test conditions: <ul style="list-style-type: none"><li>• Mobile phase: 60:40 methanol/water pre-mixed</li><li>• Flow rate: 0.5 mL/min</li><li>• Sample mix: alkylphenone mix (5.0 µL injection volume)</li></ul> |

## QSM performance specifications (Continued)

| Item                                | Specification  |
|-------------------------------------|--|
| Flow precision (continued)          | <ul style="list-style-type: none"> <li>• Column: ACQUITY UPLC BEH C<sub>18</sub>, 1.7 µm, 2.1 × 50 mm.</li> <li>• Column temperature: 35 °C ±0.3 °C</li> <li>• Detector: UV, 254 nm wavelength</li> </ul>  |
| Composition ripple (baseline noise) | <p>&lt;1.0 mAU (&lt;0.1 mAU with optional 250-µL mixer)</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• Mobile phase: A: water + 0.1%, trifluoroacetic acid; B: acetonitrile + 0.1%, trifluoroacetic acid</li> <li>• Flow rate: 0.5 mL/min</li> <li>• Gradient conditions: 1.0 to 33% B in 10 min; time average window, 10 s. Noise range 4.00 to 6.00 min</li> <li>• Column: ACQUITY UPLC BEH C<sub>18</sub>, 1.7 µm, 2.1 × 50 mm.</li> <li>• Detector: ACQUITY TUV, 214 nm wavelength, 10 Hz sampling rate</li> </ul> |
| Composition accuracy                | <p>±0.5% absolute (full scale) from 5 to 90% from 0.500 to 2.000 mL/min</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• Mobile phase: acetonitrile/water, 10:90; acetonitrile/water, 10:90, with caffeine at 12 mg/L concentration</li> <li>• Backpressure: 13,790 kPa (138 bar, 2000 psi)</li> <li>• Gradient conditions: Step gradient method</li> <li>• Detector: UV, 273 nm wavelength.</li> </ul>  |

## QSM performance specifications (Continued)

| Item                         | Specification  |
|------------------------------|--|
| Composition precision        | <p>&lt;0.15% RSD or <math>\pm 0.04</math> min SD, whichever is greater, based on six replicate injections</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• Mobile phase: 60:40 methanol/water mixed online via the QSM</li> <li>• Flow rate: 0.5 mL/min</li> <li>• Sample mix: alkylphenone mix (5.0 <math>\mu</math>L injection volume)</li> <li>• Column: ACQUITY UPLC BEH C<sub>18</sub>, 1.7 <math>\mu</math>m, 2.1 <math>\times</math> 50 mm.</li> <li>• Detector: UV, 254 nm wavelength</li> <li>• Column temperature: 35 °C <math>\pm</math>0.3 °C</li> </ul> |
| Compressibility compensation | Automatic and continuous   |
| Priming                      | Wet priming can run at flow rates up to 4 mL/min   |
| Plunger wash                 | <p>Equipped with a wash system, to flush the rear of the high pressure seal and the plunger. The interval between seal wash pump activations is:</p> <ul style="list-style-type: none"> <li>• QSM default: 5.0 min</li> <li>• bioQSM default: 0.1 min (6 s)</li> </ul>   |
| Flow ramping                 | <p>Range: 0.01 to 30.00 min to reach 2.0 mL/min</p> <p>Default: 0.45 min, to reach 2.0 mL/min at 4.44 mL/min/min</p>   |
| Vent valve                   | Used for priming the pump and automated leak testing. When the column manager switches columns, the vent valve switches to the vent position, to reduce system pressure.   |

### QSM performance specifications (Continued)

| Item              | Specification  |
|-------------------|--|
| Solvent lines     | Set of factory-installed inlet tubing assemblies. Each assembly includes a 10- $\mu$ m reservoir filter. |
| Composition range | 0.0 to 100.0% settable in 0.1% increments.   |

## Sample manager - FTN

The following table lists the performance specifications for the ACQUITY UPLC H-Class SM-FTN and H-Class bioSM-FTN.

### Sample manager - FTN performance specifications

| Item                   | Specification   |
|------------------------|---|
| Injection volume range | 0.1 to 10.0 $\mu$ L as standard. Up to 1000.0 $\mu$ L with optional extension loops (50, 100, 250, and 1000 $\mu$ L).   |
| Accuracy (aspiration)  | $\pm 0.2$ $\mu$ L, measured by fluid weight removed from vial with 10 $\mu$ L injections averaged over 20 injections using 100 $\mu$ L syringe.   |
| Linearity              | >0.999 (standard needle)<br>Test conditions: <ul style="list-style-type: none"><li>• Chromatography: isocratic</li><li>• Mobile phase: 10:90 acetonitrile/water</li><li>• Flow rate: 0.6 mL/min</li><li>• Sample mix: caffeine 0.03 mg/mL (0.2 to 10 <math>\mu</math>L injection volume)</li><li>• Column: ACQUITY UPLC BEH C<sub>18</sub>, 1.7 <math>\mu</math>m, 2.1 <math>\times</math> 50 mm.</li><li>• Detector: UV, 273 nm wavelength</li><li>• Column temperature: 40 <math>^{\circ}</math>C <math>\pm 0.3</math> <math>^{\circ}</math>C</li></ul> |


## Sample manager - FTN performance specifications (Continued)

| Item                    | Specification   |
|-------------------------|---|
| Precision               | <p>&lt;1% area RSD, 0.2 <math>\mu</math>L to 1.9 <math>\mu</math>L (0.25 to 0.50 mg/mL caffeine),<br/>         &lt;0.5% area RSD, 2.0 <math>\mu</math>L to 10.0 <math>\mu</math>L (0.03 mg/mL caffeine)<br/>         Test conditions:</p> <ul style="list-style-type: none"> <li>• Chromatography: isocratic</li> <li>• Replicates: 6</li> <li>• Mobile phase: 10:90 acetonitrile/water</li> <li>• Flow rate: 0.6 mL/min</li> <li>• Column: ACQUITY UPLC BEH C<sub>18</sub>, 1.7 <math>\mu</math>m, 2.1 <math>\times</math> 50 mm.</li> <li>• Detector: UV, 273 nm wavelength</li> <li>• Column temperature: 40 °C <math>\pm</math>0.3 °C</li> </ul>  |
| Number of sample plates | <p>Any combination (two) of the following Waters-certified plates:</p> <ul style="list-style-type: none"> <li>• 96 round-well, 350-<math>\mu</math>L</li> <li>• 96 round-well, medium, 1.0-mL</li> <li>• 96 square-well, tall, 2.0 mL</li> <li>• 384-well, 100-<math>\mu</math>L</li> <li>• 384-well, 250-<math>\mu</math>L</li> <li>• 96-well with 700-<math>\mu</math>L glass insert, with plate</li> <li>• 96-well with 1-mL extended glass insert, with plate</li> <li>• 96-well with 700-<math>\mu</math>L glass insert, insert only</li> <li>• 96-well with 1-mL extended glass insert, insert only</li> </ul> <p>For more information, see <i>Using Plates and Vials with ACQUITY UPLC and ACQUITY UPLC H-Class Systems</i>.</p> |

## Sample manager - FTN performance specifications (Continued)

| Item                                 | Specification  |
|--------------------------------------|--|
| Maximum sample capacity              | <p>768 in two, 384-well Waters-certified plates, or 96 in 2-mL vial holders.</p> <p>For more information, see <i>Using Plates and Vials with ACQUITY UPLC and ACQUITY UPLC H-Class Systems</i>.</p>  |
| Sample compartment temperature range | <p>4.0 to 40.0 °C, settable in 0.1 °C increments</p> <ul style="list-style-type: none"> <li>• Sample manager maintains setpoint temperatures from 5 to 40 °C for ambient temperatures up to 23 °C and relative humidity levels up to 80%. Sample temperatures will be within <math>\pm 3</math> °C of the setpoint.</li> <li>• At ambient temperatures above 23 °C, the sample manager can maintain an average sample temperature of 18 °C below ambient <math>\pm 3</math> °C.</li> </ul> |

Sample manager - FTN performance specifications (Continued)





## Sample manager - FTN performance specifications (Continued)

| Item   | Specification   |
|--|---|
| <b>Recommended temperature sensor locations</b><br><p>The following diagram shows the recommended temperature sensor locations on the sample tray when validating specifications.</p> <p>TS = Temperature sensor</p> |   |
| Temperature accuracy   | No more than a $\pm 0.5\text{ }^{\circ}\text{C}$ in temperature between a traceable external temperature measurement device and instrument temperature measurement device.  |
| Temperature stability  | $\pm 1.0\text{ }^{\circ}\text{C}$ (at the sensor with sample compartment door closed)   |
| Injection needle wash  | Integral, active, and programmable  |
| Minimum sample required  | <p>3 <math>\mu\text{L}</math>, residual, using 950-<math>\mu\text{L}</math> Waters Total Recovery Vials (needle placement of 0 mm) in the 48-well vial holder.</p> <p>7 <math>\mu\text{L}</math>, residual, using Maximum Recovery Vials (needle placement of 2.0 mm) in the 48-well vial holder.</p> |

## Sample manager - FTN performance specifications (Continued)

| Item                  | Specification  |
|-----------------------|--|
| Sample carryover - UV | <p data-bbox="753 253 1053 282">&lt;0.004% caffeine (UV)</p> <p data-bbox="753 291 968 321">Test conditions:</p> <ul data-bbox="761 335 1268 1269" style="list-style-type: none"> <li data-bbox="761 335 1096 364">• Solvent A: 100% water</li> <li data-bbox="761 378 1172 407">• Solvent B: 100% acetonitrile</li> <li data-bbox="761 421 1210 487">• Weak wash: water/acetonitrile, 90/10</li> <li data-bbox="761 501 1225 567">• Strong wash: water/acetonitrile, 90/10</li> <li data-bbox="761 581 1268 647">• Column: ACQUITY UPLC BEH C<sub>18</sub> 1.7 µm, 2.1 × 50 mm</li> <li data-bbox="761 661 1239 727">• Mobile phase: 90% solvent A:10% solvent B</li> <li data-bbox="761 741 1089 770">• Flow rate: 0.6 mL/min</li> <li data-bbox="761 784 1268 954">• Sample: caffeine, at 0.16 mg/mL (standard) and 4 mg/mL (challenge) in water/acetonitrile, 90/10, compared to blanks of water/acetonitrile, 90/10</li> <li data-bbox="761 968 1089 998">• Injection volume: 5 µL</li> <li data-bbox="761 1012 1160 1041">• Column temperature: 40 °C</li> <li data-bbox="761 1055 1268 1152">• Detection: UV at 273 nm, sampling rate = 20 points/s, filter time constant = normal (0.2 s)</li> <li data-bbox="761 1166 1011 1196">• Run time: 2 min</li> <li data-bbox="761 1209 1139 1269">• Data system: Empower or MassLynx software</li> </ul> <p data-bbox="753 1282 1268 1489">Basis of calculation: Any peak in the blanks following the challenge sample are compared to the known (0.005%) standard. Carryover peak areas below the standard area are within specification.</p> |

## Sample manager - FTN performance specifications (Continued)

| Item                  | Specification   |
|-----------------------|---|
| Sample carryover - MS | <p>&lt;0.005% sulphadimethoxine (MS)</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• Solvent A: water with 0.1% formic acid</li> <li>• Solvent B: acetonitrile with 0.1% formic acid</li> <li>• Weak wash: water/acetonitrile, 95/5</li> <li>• Strong wash: water/acetonitrile, 50/50</li> <li>• Mobile phase: 80% solvent A:20% solvent B</li> <li>• Flow rate: 0.3 mL/min</li> <li>• Sample: sulphadimethoxine at 5 pg/μL (standard) and 1 μg/μL (challenge) in water/acetonitrile, 90/10 +0.1% formic acid, compared to blanks of water/acetonitrile, 90/10 +0.1% formic acid</li> <li>• Injection volume: 5 μL</li> <li>• Column: ACQUITY UPLC BEH C<sub>18</sub> 1.7 μm, 2.1 × 50 mm</li> <li>• Column temperature: 40 °C</li> <li>• Sample temperature: 10 °C</li> <li>• Detection: MS SIR at 311.3 Da, 0.5 s dwell or MRM at 156.0 to 310.0</li> <li>• Ion mode: ES+</li> <li>• Run time: 5 min</li> <li>• Data system: Empower or MassLynx software</li> </ul> |

### Sample manager - FTN performance specifications (Continued)

| Item                              | Specification   |
|-----------------------------------|---|
| Sample carryover - MS (continued) | Basis of calculation: Any peak in the blanks following the challenge sample are compared to the known (0.005%) standard. Carryover peak areas below the standard area are within specification. |

## Column heater

The following table lists the performance specifications for the ACQUITY UPLC H-Class CH-A, H-Class 30-cm column heater with active pre-heater (CH-30A), H-Class bioCH-A, and H-Class bioCH-30A.

### Column heater performance specifications

| Item                                 | Specification  |
|--------------------------------------|--|
| Column capacity                      | CH-A:<br>Single column, up to 4.6-mm internal diameter (ID), to 150-mm length, with filter or guard column. Maximum column, outside diameter (OD), is 5/8-inch.<br><br>CH-30A:<br>Single column, up to 4.6-mm internal diameter (ID), to 300-mm length, with filter or guard column. Maximum column, outside diameter (OD), is 5/8-inch. |
| Column compartment temperature range | CH-A/CH-30A:<br>20 to 90 °C, in increments of 0.1 °C (control requires a setpoint of greater than ambient temperature +5 °C)   |

## Column heater performance specifications (Continued)

| Item                                     | Specification   |
|--|---|
| Column compartment temperature accuracy  | <p>CH-A/CH-30A:<br/> Tested to <math>\pm 0.5</math> °C</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• Door closed</li> <li>• No column installed</li> <li>• No flow</li> <li>• Measurement taken with traceable, external temperature measurement device</li> <li>• Measurement taken after 1 hour of thermal equilibration at set point</li> <li>• Measurement taken at column compartment sensor location</li> <li>• Tested at 35 °C, 55 °C, and 85 °C</li> </ul> |
| Column compartment temperature stability | <p>CH-A/CH-30A:<br/> Tested to <math>\pm 0.3</math> °C</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• Door closed</li> <li>• No column installed</li> <li>• No flow</li> <li>• Measurement taken with traceable, external temperature measurement device</li> <li>• Measurement taken for 1 hour of thermal equilibration at set point</li> <li>• Measurement taken at column compartment sensor location</li> <li>• Tested at 35 °C, 55 °C, and 85 °C</li> </ul>   |

### Column heater performance specifications (Continued)

| Item                 | Specification   |
|----------------------|---|
| Solvent conditioning | CH-A:<br>Active pre-heating is standard. Passive pre-heating is allowed for legacy support.<br>CH-30A:<br>Active pre-heating is standard. |

## Column manager

The following table lists the performance specifications for the ACQUITY UPLC H-Class CM-A, H-Class auxiliary column manager (CM-Aux), H-Class bioCM-A, and H-Class bioCM-Aux.

### Column manager performance specifications

| Item  | Specification   |
|---|---|
| Columns   | Up to 4.6 mm ID; up to 150-mm columns. Up to 30-mm column with guard column or filter. ACQUITY UPLC 30-cm column heater/cooler supports columns with dimensions of up to 7.8 mm diameter.                       |
| Number of columns                               | CM-A:<br>Up to two 150-mm columns (with filters or guard column) or up to four 50-mm columns (with no guard column or filter).<br>CM-Aux:<br>Up to two 150-mm columns (with filters or guard column).           |
| Column compartment temperature range (settable) | 4 to 90 °C, in increments of 0.1 °C. Troughs are independently settable.<br>Derating: The minimum achievable column compartment temperature set point must not be greater than 25 °C below ambient temperature. |

## Column manager performance specifications (Continued)

| Item   | Specification   |
|--|---|
| Time to temperature, from steady state, after door is open for 30 seconds. | <p>12 minutes maximum</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• No column installed</li> <li>• No flow</li> <li>• Measurement taken with internal temperature sensor</li> <li>• Measurement taken after 1 hour of thermal equilibration at set point</li> <li>• Door is opened for 30 seconds</li> <li>• Tested at 35 °C, 55 °C, and 85 °C</li> </ul>  |
| Column compartment temperature accuracy                                    | <p>Tested to <math>\pm 0.5</math> °C</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• Door closed</li> <li>• No column installed</li> <li>• No flow</li> <li>• Measurement taken with traceable, external temperature measurement device</li> <li>• Measurement taken after 1 hour of thermal equilibration at set point</li> <li>• Measurement taken at column compartment sensor location</li> <li>• Tested at 35 °C, 55 °C, and 90 °C</li> </ul> |

## Column manager performance specifications (Continued)

| Item  | Specification   |
|---|---|
| Column compartment temperature precision          | <p>Tested to <math>\pm 0.1</math> °C</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• Door closed</li> <li>• No column installed</li> <li>• No flow</li> <li>• Measurement taken with traceable, external temperature measurement device</li> <li>• Measurement taken at column compartment sensor location</li> <li>• Temperature is ramped from ambient to 90 °C</li> <li>• Measurement taken after 1 hour of thermal equilibration</li> <li>• Temperature is returned to ambient</li> <li>• Test is repeated four additional cycles</li> </ul> |
| Column compartment temporal temperature stability | <p>Tested to <math>\pm 0.3</math> °C</p> <p>Test conditions:</p> <ul style="list-style-type: none"> <li>• Door closed</li> <li>• No column installed</li> <li>• No flow</li> <li>• Measurement taken with traceable, external temperature measurement device</li> <li>• Measurement taken for 1 hour after thermal equilibration at set point</li> <li>• Measurement taken at column compartment sensor location</li> <li>• Tested at 35 °C, 60 °C, and 90 °C</li> </ul>  |
| Ambient temperature stability                     | Within 2.0 °C/60 minute maximum   |
| Pre-heater temperature (not user-settable)        | Defined by set point of the column trough.  |



## Wetted materials of construction

### Quaternary solvent manager

The following table lists the wetted materials of construction for the ACQUITY UPLC H-Class QSM and H-Class bioQSM.

#### QSM wetted materials of construction

| Description      | Specification   |
|------------------|---|
| Wetted materials | ACQUITY UPLC H-Class QSM:<br>316 L stainless steel, DLC, fluoroelastomer, fluoropolymer, Nitronic 60, PEEK™ and PEEK blend, PPS and PPS blend, ruby, sapphire, titanium alloy, UHMWPE blend, zirconia<br>ACQUITY UPLC H-Class bioQSM:<br>DLC, fluoroelastomer, fluoropolymer, Inconel 600, MP35N, PEEK™ and PEEK blend, PPS and PPS blend, ruby, sapphire, titanium, titanium alloy, UHMWPE blend, zirconia |

### Sample manager - FTN

The following table lists the wetting materials of construction for the sample ACQUITY UPLC H-Class SM-FTN and H-Class bioSM-FTN.

#### Sample manager - FTN wetted materials of construction

| Description      | Specification  |
|------------------|--|
| Wetted materials | ACQUITY UPLC H-Class SM-FTN:<br>316 stainless steel, 316 stainless steel with DLC, 450-G, borosilicate, EPDM, fluoropolymer, fluoropolymer-coated Neoprene, gold-plated stainless steel, PEEK, PEEK blend, polyethylene, polyimide, PPS<br>ACQUITY UPLC H-Class bioSM-FTN:<br>MP35N, PEEK, titanium alloy with DLC |

## Column heater

This table lists the wetted materials of construction for the ACQUITY UPLC H-Class CH-A, H-Class CH-30A, H-Class bioCH-A, and H-Class bioCH-30A.

### Column heater wetted materials of construction

| Description      | Specification  |
|------------------|--|
| Wetted materials | ACQUITY UPLC H-Class CH-A and CH-30A:<br>316 stainless steel<br>ACQUITY UPLC H-Class bioCH-A and bioCH-30A:<br>MP35N |

## Column manager

This table lists the wetted materials of construction for the ACQUITY UPLC H-Class CM-A, H-Class CM-Aux, H-Class bioCM-A, and H-Class bioCM-Aux.

### Column manager wetted materials of construction

| Description      | Specification  |
|------------------|--|
| Wetted materials | ACQUITY UPLC H-Class CM-A and CM-Aux:<br>316 stainless steel<br>ACQUITY UPLC H-Class bioCM-A and bioCM-Aux:<br>MP35N, titanium |